

ACUTECH BLIND QUALITY CONTROLS

Acutech specimens are designed to comply with the blind testing requirements and objectives of the Mandatory Guidelines for Federal Workplace Drug Testing Programs. Protox Services can supply your facility certified samples and assist you and your staff in complying with blind sample regulations.

Acutech specimens are manufactured in a human urine matrix and are certified by SAMHSA Certified Laboratories for drug content. These performance samples are packaged in 60 ml plastic bottles and are shipped as a liquid, ready-to-use product. Acutech performance samples are processed as a normal patient collection. Select the blind sample to be tested from the refrigerator, remove the cap and inner seal, and pour Acutech into your laboratory's transport container(s). The labels on all Acutech bottles are removable and attach to the MRO Copy of the Chain of Custody.

Acutech specimens are sold in any quantities or combinations including certified negative, positive, dilute, substituted and adulterated samples for any regulated workplace drug testing program you may be enrolled in.